

KJMR-II Blood Roller Mixer

Description:

The product is used for mixing and blending of the blood, powder reagents, total blood control materials in the anticoagulation tube or bottle. It is proved of being more accurate and effective than the manual mixing by the Shanghai Clinical Inspection Center as well as the Big Hospital all over that China. According to different requirements, it can be rolling mixing, anticlockwise rolling and speed adjustable. It is an ideal instrument for clinica I Lab.



Speed:

Vibrating: 10~80R

PMRolling: 10~80

RPM

Power consumption: 12W

KJMR-IIA Blood Roller Mixer

Description:

KJMR-IIA is the upgrade model of KJMR-II. Digitally controlling on the speed and time is applied to this new model on the basis of the old model function which is clearly displayed by LED. It would be more convenient for the operator and improve the accuracy of the experiment.

Spec: Speed swing: 20-100times/minRolling: 20-100r/min

Swing amplitude: 22 ± 1 mm

Timing range: 0-999min

Running mode: Timing/continuous

Power rating: 10W

Voltage: AC100-240V(50/60Hz)

N.W.: 3.6Kg

External dimension: 394*266*98mm

